

## **SSOCTA artefact removal and image enhancement**

Carl Glittenberg, MD

Karl Landsteiner Institute for Retinal Research and Imaging,  
and Topcon Europe Medical B.V.

This presentation explores the source of artefacts in OCT-A and how to remove them. White line removal, black band removal, speckle noise removal, horizontal flicker noise removal, and projection artifact removal will be discussed. Additionally, methods for enhancing the quality and information content of OCT-A imagery will be discussed as will the benefits of 3D visualization. The presentation shows that when the data is corrected and enhanced in the right way, considerably more information can be gathered than initially apparent.